

Chemical Casserole

Purpose

To show different types of reactions (precipitation reaction and complex ion formation)

Materials

6 beakers	17M NH_4OH
1M AgNO_3	1M KBr
0.1M NaOH	1M $\text{Na}_2\text{S}_2\text{O}_3$
1M NaCl	1M KI

Procedure

- 1) Mix 5mL 1M AgNO_3 and 20mL H_2O
- 2) Fill each beaker with one of the following solution of the given volume:
 - a. 10mL of 1M KI
 - b. 10mL of 0.1M NaOH
 - c. 10mL of 1M NaCl
 - d. 15mL of 17m NH_4OH
 - e. 10mL of 1M KBr
 - f. 20mL of 1M $\text{Na}_2\text{S}_2\text{O}_3$
- 3) Add just enough of the AgNO_3 solution into each beaker as to allow reactions to occur. Solutions should be mixed in the following order:
 - a. 10mL 0.1M $\text{NaOH} \rightarrow \text{AgOH}$ (brown)
 - b. 10mL 1M $\text{NaCl} \rightarrow \text{AgCl}$ (white)
 - c. 15mL 17m $\text{NH}_4\text{OH} \rightarrow \text{Ag}(\text{NH}_3)_2^+$ (clear)
 - d. 10mL 1M $\text{KBr} \rightarrow \text{AgBr}$ (off-white)
 - e. 20mL 1M $\text{Na}_2\text{S}_2\text{O}_3 \rightarrow \text{Ag}(\text{S}_2\text{O}_3)_2^{3-}$ (clear)
 - f. 10mL 1M $\text{KI} \rightarrow \text{AgI}$ (yellow)

Additional Information

- 1) Possible reactions can involve cyanide and sulfide ions however CN^- is too hazardous and S^{2-} is too smelly

- 2) At any point, the demonstration may be stopped and equilibrium constants used to explain an observation or to predict an upcoming change
- 3) Just enough solution can be added at each point to create or dissolve the precipitate

Disposal

Solutions should be placed in properly labeled waste containers with UI# 100963.